

Functional Anaesthetic Discography™ Procedure Pain Score Card

Patient Data: _____

Date: _____

Level 1 _____

Type of Anesthetic: _____ Volume of Anesthetic: _____

NPRS Score	Baseline	During Injection	After Injection	Other
Time				
Pain Eliciting Activity				
Forward Flexion				
Other Postures (optional)				

Range of Motion – Forward Flexion

Mark the location of hands during forward flexion.

Level 2 _____

Type of Anesthetic: _____ Volume of Anesthetic: _____

NPRS Score	Baseline	During Injection	After Injection	Other
Time				
Pain Eliciting Activity				
Forward Flexion				
Other Postures (optional)				

Range of Motion – Forward Flexion

Mark the location of hands during forward flexion.

Level 3 _____

Type of Anesthetic: _____ Volume of Anesthetic: _____

NPRS Score	Baseline	During Injection	After Injection	Other
Time				
Pain Eliciting Activity				
Forward Flexion				
Other Postures (optional)				

Range of Motion – Forward Flexion


Mark the location of hands during forward flexion.

Notes: _____

1 Have patient elicit typical back pain. Record patient's pain score(s) and range of motion



2 Inject 1 cc of anesthetic into the targeted disc




3 Have patient reassume the same pain-eliciting posture or activity. Record patient's pain score(s) and range of motion



Notes:

4 Compare pre- and post-injection pain scores and range of motion to identify source(s) of discogenic pain

Level 1					L4-S				Range of Motion - Forward Flexion			
Type of Anesthetic	4% Lidocaine			Volume of Anesthetic:	1 cc				Baseline	After Injection	Baseline	After Injection
NPRS Score	Baseline	During Injection	After Injection	Other								
Time	1:05	1:08	1:18					<input checked="" type="checkbox"/>				
Pain Eliciting Activity								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<i>Sitting</i>	8	10	2					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Forward Flexion	7		2					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Other Postures (optional)								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<i>Sit to Stand</i>	8		2					<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



5 Connect a syringe to the stopcock assembly. Open the stopcock and deflate the balloon by pulling a vacuum with the syringe, holding for 5 seconds. Leaving the stopcock open, slowly remove the catheter(s) from the patient

MEDTRONIC
Spinal and Biologics
1221 Crossman Avenue
Sunnyvale, CA 94089
USA
877.459.7466
Tel: 408.548.6500
Fax: 408.548.6501
www.Discyphor.com

As with most interventional procedures, the Functional Anaesthetic Discography™ (F.A.D.™) Procedure has associated risks, including serious complications. For complete information regarding indications for use, contraindications, warnings, precautions, adverse events and methods of use, please reference the devices' Instructions for Use.

© 2008 Medtronic Spine LLC. All Rights Reserved. 16000756-02